# Ena Flow Sealant

## (IT) ITALIANO

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# Ena Flow Sealant

## (EN) ENGLISH

Ena Flow Sealant is a light-cure, pit and fissure sealant.

### Composition:

Monomer Matrix: Diurethandimethacrylate; Tetramethylene dimethacrylate  
Total content of fillers: 33 % weight, inorganic fillers (0,005-3,0 µm)

### Indications

Sealing the enamel pits and fissures of teeth of different form to aid in the prevention of caries.

### Contraindications

In case of known allergy to some of the components do not use it

### Precautions

SEALANT: Contains Tetramethylene dimethacrylate, Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide

Warning: May cause an allergic skin reaction. Wear protective gloves/protective clothing/eye protection/face protection. If skin irritation or rash occurs: Get medical advice/attention.

## INSTRUCTIONS FOR USE

**Note:** Sealant is light-sensitive.

**Teeth selection:** Select teeth sufficiently erupted so that a dry field can be maintained. The morphology of the pits and fissures should be deep.

**Preparation:** Clean with fluoride-free prophylaxis paste. We suggest using a rubber dam to isolate the area.

**Etching:** Etching according to manufacturer's instruction.

**Application of Sealant:** Take Ena Flow Sealant out from syringe using the application tip and a brush (Enamel plus M brush) and slowly introduce it into the pits and fissures. Do not let sealant flow beyond the etched surfaces. Stirring the sealant with the syringe-tip during or after placement will help to eliminate possible bubbles, and to enhance the flow into the pit and fissures. Cure the sealant by exposing it to light from a light curing unit (CLEDPPLUS Micerium), or other curing unit of comparable intensity ( $\geq 1200$  mW) for 20 seconds, from all sides of the tooth; keep the light-curing tip as close as possible to the restoration. When set, the sealant forms a hard, opaque film, off-white in colour with a slight surface inhibition. Inspect sealant to make sure that coverage is complete and there are no gaps. If surface has been contaminated, additional sealant should be added. If contamination has occurred re-etch, rinse, and dry prior to placing more sealant.

Wipe the sealant with a cotton swab to remove the thin sticky film on the surface. Check occlusion and adjust as required.

**Note: instructions for Flow syringe.** Apply the application tip on the syringe after unscrewing the cap. Being composite flowable, when you push the piston you will activate a thrust and the material will start and continue to come out. To stop the flow it is enough to pull the piston back only 1 mm. Careful: avoid pulling piston excessively, otherwise air can come into the syringe and air bubbles will enter in the following emissions of material. A minimum opposite movement is enough, the piston will return in position elastically, avoiding air bubbles. To avoid the excessive emission of material we suggest to hold the tips of the syringe direct upwards till next application on the same patient. We also suggest starting pushing the piston in this position, so if there would be air in the syringe, bubbles will come out before the material. At the end of the restoration, remove the application tip and place the cap on the syringe again screwing it.

### Disinfection / Protection from cross-contamination

Place the syringe with the application tip attached into a suitably shaped barrier sheath; pierce the end of the sheath, exposing the cannula for use. The use of a barrier sheath facilitates cleaning and disinfection of the syringe between patients.

After use, hold the tip through the sheath; unscrew and remove the tip along with the sheath. Dispose of the used tip and the sheath through the appropriate waste stream. Place the cap on the syringe again screwing it.

Disinfection - After removing the application tip and the sheath, disinfect the syringe using an intermediate-level disinfection process (liquid contact) as recommended by the Centre for Disease Control and endorsed by the American Dental Association. Guidelines for Infection Control in Dental Health-Care Settings - 2003 (Vol.52; No. RR-17), Centre for Disease Control and Prevention.

### Notes for light curing:

Use a light curing unit with an emission spectrum of 350-500 nm to cure this material. As the required physical properties can only be achieved if the unit works correctly, its luminous intensity must be checked regularly as described by the light curing unit manufacturer.

## USE AND STORAGE

Do not store above 25°C (77°F) and below 3°C (38°F) .

Do not use the product after the expiration date (see label on syringe).

Due to hygienic reasons flow application tips should be used only once.

Use the material at room temperature.

Medical device, for dental use only: keep away from children.

After use, close container with cap and keep it closed. Avoid direct exposure to sunlight.

If the material is not completely cured, it may discolour, mechanical properties deteriorate and pulpal inflammation can occur.

See the Material Safety Data Sheet for disposal information.

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**MICERIUM S.p.A.**