# Temp Silicon Brush

TPENMC  
TPENMS

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### **Sterilizable silicon brushes for composite application**

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## **Sterilizable silicon brushes for composite application**

- They do not stick on composite. They do not leave porosity and can be handled as normal brushes used by the lab.
- They have a resilience similar to uncured composite: this avoids air inclusion, and allows each layer of material to be compacted and smoothed.
- Remove excess soil under running water. Disassembly the rubber tip from the metal handle. Automated Cleaning: Ultrasonic Cleaner. Manual Cleaning: Detergent (alcohol), disposable tissue, running water. Sterilization in autoclave. Rapid cycle: Metal handle: 3 min. at 134°C; Rubber tips: 15 min. at 121°C.
- They can be used also with other materials: acrylics, monomers, sealers, lubricants, cianoacrylates, separators, gels, glazures, stains.

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### **INTRODUCTORY KIT TEMP BRUSHES**

5 assorted silicon tips for the modellation of composite / 1 metal handle straight + 90°

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### **INTRODUCTORY KIT TEMP BRUSHES**

5 assorted silicon tips for the modellation of composite / 1 metal handle curved + 45°

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Silicon tip short spatula shape for modelling occlusal surfaces

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Silicon tip long spatula shape for blending layers and shaping surfaces

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Silicon tip round-flat shape for modelling cervical areas

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Silicon tip cone shape point ideal for occlusal shaping

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Silicon tip with three sides for stains and characterizations

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