TENDERFLASK

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Spray the mould with Tender silicone separator and reposition the wax-up.

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TENDERFLASK

Pressing system for light curing composite

The Tender Flask System allows the reproduction in light curing resin of crown & bridge and aesthetic sections of partials, through the pressing technique in a flask. The technique is very simple and offers several advantages; operative times are considerably reduced and waiting times are limited only to the silicone and composite setting time. This system requires two different silicones:

- TEMP-SILIC PUTTY, 90 shore putty hardness silicone, for the base
- TEMP-SILIC CLEAR, 50 shore transparent silicone in self-mixing cartridge, for the mould

The wax-up is perfectly reproduced with all details, by respecting shape and function, thanks to the transparent silicone mould. No special unit is required

INSTRUCTIONS FOR USE

After casting finish and prepare the metal structure.

Proceed with the wax-up of the aesthetic elements that should be reproduced in light curing composite, respecting anatomy and function. Spray Tender Silicone separator in the flask.

Fill up the flask base and the inner part of the metal structure with the 90 shore silicone TEMP-SILIC PUTTY after mixing paste A + B (mixing ratio 5:5 +5 retarding drop; working time 1-2 min.), and place the framework in the flask base. When the silicone is set (after 2-3 min.), remove the wax-up, eliminate debris and silicone excess up to the margin area. Spray the mould with Tender silicone separator and reposition the wax-up. Then, fill up the cover of the flask with TEMP-SILIC CLEAR silicone (working time: 5 min.), and close the flask. When the silicone is set (after about 15-20 min. on the bench, or 10 min. under pressure at 2 atm dry and cold setting, max. 40°C), open the flask and remove the excess of silicone. Remove the wax-up and prepare the metal by eliminating the wax, sandblasting, and cleaning the metal as suggested by the manufacturer. Apply metal primer and opaque following the instructions of the manufacturer. Prepare some venting channels for the resin in the hard silicone of the base with a ball bur. Spray some separator on the mould, reposition the framework and apply the dentine body in the transparent mould***.

Close the Tender Flask and leave it to sit for about 8 min. in a curing chamber without pressure (dark still phase). Then insert the flask in the light curing unit increasing the suggested curing time by 4 minutes. After curing, open the flask, remove the framework, and create the enamel cuts with fine grain burs. Clean and apply a bonding (like Temp bonding fluid – ref. TEMPF) to the surface. Apply opalescence and characterization masses and light cure them. Then apply the enamel body in the mould and close the flask. Repeat dark still phase and cure again as indicated above.

Open the flask, finish and polish.

*** With Enamel plus Composite system, before pressing of the dentine, stratify the Tender masses, for strength and aesthetic reasons.

Micerium S.p.A. Via Marconi, 83 - 16030 Avegno (GE)
Tel. 0185 7887 870 - Fax 0185 7887 970
www.micerium.it - e-mail:micerium@micerium.it