ENA  *post*

Perno endocanalare in fibra di vetro
Glass Fiber Post
Glasfiber-Wurzelstift

estefica  [®]
GRUPPO
MICERIUM

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Dentista

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(EN) Glass fiber post - ENA POST

Perfect aesthetic and functional integration

Benefits and features

- Perfect cohesion fibers-composite: fibers are inserted in the same resin matrix of the microhybrid composite for restoration
- Aesthetic: dentine colour and fluorescent
- Functional: same elasticity of dentine
- Optimal X-ray diagnosis: radiopaque

Intended purpose

Reinforcement of restoration of endodontically treated teeth

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
*The coloured rings can be used to measure the depth of the canal, positioning it where the post must be cut. **Note: the rings must be removed before cementation.***

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Carbon, Standard Glass fiber post and Ena Post glass fiber post (right picture under UV light).

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Ena Post with 2% conicity (right picture under UV light).

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Ena Post with 10% conicity (right picture under UV light).

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Ena Post is a post made of glass fibers and resin used to reinforce restorations or core build-up in devitalized teeth. It is available in 2 conicities (2% and 10%) and 4 sizes each. It is used after preparing the root canal with a colour-coded drill calibrated with the post.

Patient target group and medical condition

Medical device intended for patients who have been treated with devitalization of a tooth, and now need a dental restoration: Children 6-18 years (permanent teeth only), adults 19-64 years, elderly 65- above, of any sex and condition.

Indications

Ena Post is used for:

- Reinforcement of aesthetic restorations in devitalized teeth
- Reinforcement of core build-up in devitalized teeth

Note: Any serious incident relating to the device must be reported to the manufacturer (Micerium S.p.A.) and to the competent authority of the Member State in which the user and/or patient is established.

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REPROCESSING INSTRUCTIONS FOR REUSABLE MEDICAL DEVICES (ISO EN 17664)

Device(s): calibrated drills for Ena Post.

WARNINGS

Instruments are supplied not sterile and should be sterilized before the first use and after each use. Appropriate handling precautions should be taken during use and reprocessing. Gloves, eye protection and a mask should be worn.

Limitations on reprocessing

Repeated reprocessing has minimal effect on these instruments; the working part is deteriorating mainly with use.

INSTRUCTIONS

Point of use: Remove excess soil under running water.

Container: No particular requirements. They should be protected from damage to the working part: burs holder is recommended

Transportation: It is recommended that instruments are reprocessed as soon as possible after using them. If transported wet there is an increased chance of staining or corrosion. Prolonged storage in disinfectant solutions may result in corrosion and should be avoided, however, these solutions must be validated for use with rotating instruments.

Cleaning:

- **automated cleaning:** Ultrasonic Cleaner (e.g. Biosonic) with general purpose liquid type UC30 or other liquids validated for cleaning of dental rotary instruments (manufacturers' instructions must be followed). Specific disinfectant solutions for ultrasonic can also be used only if validated for use with rotating instruments, specifically, for brushes and rubber points.

1. Load burs directly in a beaker or in a burs holder.

2. Run cycle, minimum 15 minutes.

3. When unloading check instruments and tips for complete removal of visible soil. If necessary repeat cycle or use manual cleaning.

- **manual cleaning:** 1. Rinse soil in excess from instrument under running cold water. 2. keep the instruments immersed in a sink reserved for the purpose in which case a disinfectant agent validated for processing of dental rotary instruments must be used and the disinfectant manufacturers' instructions must be followed and use wire brushes to remove any contamination. Care should be taken to avoid spreading contaminants by spraying or splashing during the brushing process 3. Rinse under clean running water or, even better, with ionized water.

Disinfection: Disinfectant solution may be used in accordance with label instructions. Only disinfectant agent validated for processing of dental rotary instruments must be used.

Drying: When drying is achieved as part of a washer disinfectant cycle do not exceed 140°C.

Maintenance: No particular requirements

Inspection and Function - Testing: After clearing, inspect all instruments, with the aid of magnification if necessary, to ensure that all contamination has been removed. Repeat the cleaning process if necessary. Visually inspect for damage and eliminate broken or deteriorated instruments.

Packaging: single: a standard packaging material may be used. In sets: Instruments may be loaded into dedicated bur holder, wrap bur holder using appropriate method.

Sterilization: Vacuum autoclave (ISO17665). Cycle: 3 minutes minimum (18 minutes for diamond and carbide burs and Ena Post drills), 134°C maximum, 2,1 bar (follow autoclave manufacturers instructions).

Storage: No specific limitations related to storing duration. Storage should be in dry, clean conditions and at ambient temperature containers or pouches must be dry before opening to avoid recontamination of the contents with water.

Additional Information: When sterilizing multiple instruments with one autoclave cycle ensure that the sterilizer's maximum load, suggested by the manufacturer, is not exceeded.

Manufacturer contact: See brochure for telephone and address of local representative or telephone (+39) 0185 7887880

The instructions provided above have been validated by the medical device manufacturer as being CAPABLE of preparing a medical device for re-use. It remains the responsibility of the processor to ensure that the processing as actually performed using equipment, materials and personnel in the processing facility achieve the desired result. This requires validation and routine monitoring of the process. Likewise any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

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Extracted tooth endodontically treated. UV light shows again a good fluorescence.

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Tooth after 1 hour treatment of 35% hydrogen peroxide (Ena White). UV light shows lower fluorescence.

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Composite restoration with Ena Post Fiber Post. UV light shows a better fluorescence than before restoration.

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Adhesive relation: dentine - enamel dentine adhesive - composite cement - post.

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Section of restored tooth using Ena Post conicity 10% and Enamel plus (Ena HRI) composite (picture under UV light in the middle, transilluminated section on the right).

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Section of restored tooth using Ena Post conicity 2% and Enamel plus (Ena HRI) composite (transilluminated section on the right).

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Fig. 50.
 Rx finale.
Final X-Ray.
Abschluss-Röntgenbild.

(EN) Instructions for Use

A. Complete the endodontic therapy, and prepare the tooth for the restoration by evaluating if a post is necessary, based on residual quantity of intact tooth. The loss of one or both interproximal walls in anterior teeth, of one or both mesial and distal cuspid in posteriors, suggest you to use a post, in order to better stabilize the restoration.

B. Make an X-ray to determine the appropriate diameter and depth of post space preparation. The post diameter should be slightly lower than the endodontic preparation. The post should be placed at same radicular depth as the height of dentinal body (for composite aesthetic restoration) or of post and core. **Use of rubber dam is recommend.**

C. Remove gutta-percha to preplanned depth with a Gates-Glidden drill or Peeso reamer. Radiographic verification is recommended.

D. Select the Ena Post Drill diameter which corresponds to the last Gates-Glidden drill to begin preparing the post space. Use Drills with a slow-speed contra-angle with water spray in order to avoid tissues overheating. The drill must be kept in continuous clockwise rotation until it has been completely removed from the tooth. This will minimize the risk of the drill seizing in the post space. Remove any tooth debris from canal by irrigating post space with water spray. Sequentially step up to the next larger **Ena Post Drill** until the preplanned diameter and depth are achieved.

E. Select the **Ena Post** that corresponds to the last drill used to prepare the post space and insert post into post space. For a better fitting shorten the post as necessary from either the apical or occlusal end (as clinical judgement dictates) with a Diamond separating disc, wetting the post with water to avoid fibers overheating. The coloured rings can be used to measure the depth of the canal, positioning it where the post must be cut. Note: the rings must be removed before cementation.

F. Remove the post and gently clean it with an alcohol wipe and apply on the surface a mixture of **Ena Bond** with **Ena Bond Catalyst** (one drop+one drop).

G. Etch the cavity (coronal and radicular preparation) with Ena Etch 37% phosphoric acid for 2 minutes. For ideal adhesion, before etching, a microblasting on cavity surfaces is recommended in order to clean and eliminate endodontic material debris.

H. Wash accurately the canal with a syringe to remove completely the acid. Suck water and dry the canal with paper point; in order to maintain wet the dentine and avoid collagen collapse do not dry with air.

I. Apply in the cavity and in the canal the mixture of **Ena Bond** and **Ena Bond Catalyst**.

Ricostruzione di perno moncone con Ena Post (conicità 10%)
Restoration of Post and core with Ena Post (conicity 10%)
Restauration durch Stift-/Stumpf Aufbau mit Ena Post (10% Konizität)

Fig. 51.
Moncone protesico gravemente compromesso
Severly compromised prothetic core
Äusserst geschwächter prothetischer Stumpf

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The adhesive should be wiped on the surface with a disposable microbrush or with a paper point; dry with air to eliminate water and solvent residuals. Insert the post to check the canal and better push adhesive in dentine tubules.

J. Apply a disposable automixing tip and extract **Ena Cem** dual composite cement from the syringe: the two components will be automatically mixed. Place the cement directly into the canal filling it from the bottom of the cavity to the surface using a disposable endodontic applicator (Hawe-Neos); do not use any lentulo or rotating instrument. Apply some cement on the post surface and insert slowly the post to full depth, allowing excess cement to vent. Find the right position and, maintaining the post in position, wipe away any excess cement.

Note: Ena Cem working time 4 min. Complete setting time 7-9 min.

K. Light cure for 60 seconds and proceed to the restoration using **Enamel plus** (Ena HRI) composite, both for aesthetic direct restoration and post and core, that could be immediately prepared once the restorative procedure has been completed.

Warnings. Posts are single use: re-use compromises the structural integrity. Posts are supplied non-sterile. Posts can be sterilized **only** by cold sterilant.

Hazard statement: Contains tetramethylene dimethacrylate. May cause an allergic skin reaction.

Disposal: parts and accessories in direct contact with patient's mouth must be sterilised before disposal or disposed of as special waste. Disposal of the medical device must be carried out in accordance with local regulations. Contaminated packaging can be disposed of, after cleaning, in the separate collection of rubbish in accordance with the identification symbols, if applicable (97/129 EC).

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MICRIUM S.p.A. Via G. Marconi, 83 -16036 Avegno (GE) Italy