

Ena White 2.0



INNOVATIVE
BRUSHING
TECHNIQUE
WHICH
ACCELERATES
BLEACHING



ENA® WHITE 2.0

The fastest professional home bleaching with in office pre-treatment

is a system for cosmetic bleaching that uses a 6%-hydrogen-peroxide-based bleaching gel with the special XS151™ accelerator, which activates by brushing; this action exponentially increases the absorption rate of peroxide if compared to traditional bleaching.

TECHNOLOGICAL & OPERATING ADVANTAGES

- Easier to use thanks to the toothbrush with dispenser
- Greater comfort for the patient, two minutes correspond to one night with conventional, annoying bleaching trays. No ingestion of hydrogen peroxide.
- Less risk of irritation by reducing contact time: 2/3 min. per day instead of 6/8 hours to achieve the same bleaching result with bleaching trays.
- A conventional 20-day treatment provides approx. 120/160 hours contact time compared to 40/60 min. with Ena®White 2.0.

ENA WHITE 2.0 Pre-filled



ECW01



ECW01B



Ref	Description	Q.ty
ECW01*	ENA WHITE 2.0 Professional home bleaching system Toothbrush filled with 6ml H ₂ O ₂ gel	1 pc
ECW01B*	ENA WHITE 2.0 Professional home bleaching system Toothbrush filled with 6ml H ₂ O ₂ gel Enamel plus toothpaste	1 pc 1 pc
ECW02	ENA WHITE 2.0 Professional home bleaching system Toothbrush filled with 6ml H ₂ O ₂ gel	2 pcs

*minimum order: multiple of 24 pcs packed in dispenser of 12 tubes
**toothbrushes delivered in separate box for easier fridge stock (see below picture)

PHASE 1 • IN-OFFICE



PHASE 2 • AT HOME



2 PHASES

1st In-office, where mechanical brushing is performed;
2nd At home with the special Ena®White 2.0 toothbrush with dispenser filled with bleaching gel (enough for a treatment of about 20 days with 2 daily applications).



In-Office operating protocol

- Tartar removal
- Baking soda application with Smart Prophy Pink brush, mounted on contra-angle handpiece (200/250 rpm)
- Application of Ena®White 2.0 gel: leave on for 10 to 15 min.
- Brush the bleaching gel with dry Smart Prophy Pink brush, mounted on contra-angle (max 300/400 rpm) for about 30 to 40 seconds on each tooth.
- Filling toothbrush dispenser with Ena®White 2.0 gel from syringe
- After 20 days the toothbrush dispenser can be filled again with bleaching gel

Home Bleaching procedure

- Daily oral hygiene
- Spill a small amount of gel on the toothbrush and brush the teeth for 30/45 sec.
- Rinse and repeat the treatment a second time
- It is possible to continue the bleaching treatment after 20 days refilling the toothbrush with the same syringe during the follow-up visit.
- Additional bleaching with further refills are possible, with minimum intervals of 3 months

ENA WHITE 2.0
To be filled



ECW01BS



ECW6ES

Ref	Description	Q.ty
ECW01BS	ENA WHITE 2.0 Professional home bleaching system for a treatment of about 40 days	1 pc
	Toothbrush with dispenser (empty)	11 ml
	Refill H ₂ O ₂ gel	1 pc
	Enamel plus toothpaste 75 ml	

Ref	Description	Q.ty
ECW6ES	ENA WHITE 2.0 Professional home bleaching system for a treatment of about 40 days	6 pcs
	Toothbrush with dispenser (empty)	
	Refill H ₂ O ₂ gel (economy package without tubes)	6 pcs x 11 ml

Other packages available on request