

BIO-ACTIVATORS DTM • SECOND & THIRD CLASS



MODEL SC

Appliance for the treatment of Angle's second skeletal classes. Corrects second skeletal classes by repositioning the dental arches in the key of first class resulting in reconditioning of the temporomandibular joints (TMJ).



CLASS II INDICATIONS

- used in mixed deciduous and permanent dentition
- corrects skeletal class IIs by repositioning the dental arches in a class I key
- corrects midlines and balances TMJs
- corrects dysfunctional vicious habits
- results in pre-alignment of the teeth



MODEL TC

Appliance for the treatment of Angle's third skeletal classes. Corrects third skeletal classes by repositioning the dental arches in the key of first class resulting in reconditioning of the temporomandibular joints (TMJ).



CLASS III INDICATIONS

- used in mixed deciduous and permanent dentition
- corrects skeletal class IIIs by repositioning the dental arches in a class I key
- corrects midlines and balances TMJs
- corrects dysfunctional vicious habits
- results in pre-alignment of the teeth

CONTRAINDICATIONS

- dental elements with accentuated mobility of periodontal character
- overt gingivitis

MEASURES & CODES

MODEL	1/45 mm	2/50 mm	3/55 mm	4/60 mm	5/65 mm
SC	BIO-BSC1DTM	BIO-BSC2DTM	BIO-BSC3DTM	BIO-BSC4DTM	BIO-BSC5DTM
TC	BIO-BTC1DTM	BIO-BTC2DTM	BIO-BTC3DTM	BIO-BTC4DTM	BIO-BTC5DTM