



# COMPOSITE & OPAQUE FOR GUM LIGHT-CURING INDIVIDUAL TRAY & RESIN STAINS



Partial



Application of pink opaque



Opaque composite



Combinated prostheses in Tender composite and resin

## TENDER PINK COMPOSITE

Light curing composite pink for gingival flange  
Pink opaque for partial



TENDERPKIT

Ref	Description	Q.ty
<b>TENDERPKIT</b>	<b>TENDER PINK KIT</b> Light curing composite Pink Pink Orange, Pink Light, Pink Dark, Pink Transparent Pink paste opaque Enamel plus Stain Flow: RED, BLUE, WHITE	4x2,5 g  1 g 3x2 g
<b>ETPK</b>	<b>TENDER PINK KIT</b> Light curing composite Pink Pink Orange, Pink Light, Pink Dark, Pink Transparent	4x2,5 g
<b>ETP</b>	Tender Pink Orange	5 g
<b>ETPL</b>	Tender Pink Light	5 g
<b>ETPD</b>	Tender Pink Dark	5 g
<b>ETPT</b>	Tender Pink Transparent	5 g
<b>ETOP</b>	Pink paste opaque light curing for covering partial	3 g
<b>ETFP</b>	Tender Pink Orange Flow	2 g

## TEMP SILIC SILICONE

Special silicones for stent, ideal also for transparent flask



Ref	Description	Q.ty
<b>TEMPSMK</b>	<b>SILICONE KIT FOR TENDER FLASK</b> Silicone for stent/flask 90 shore A+B+retarding Silicone for stent clear 50 shore cartridge Silicone-flask separator	300 g 50 ml 60 ml
<b>TEMPSMR3</b>	Silicone for stent and flask 90 shore Paste A grey, Paste B brown + retarding	1,5+1,5 Kg
<b>TEMPSM10</b>	Silicone retarding for flasking technique	10 ml
<b>TEMPSMT</b>	Silicone for stent clear 50 shore cartridge working 1'30"-2'30", hardening 15-20 min.	50 ml
<b>TEMPSMT5</b>	Silicone for stent clear 50 shore cartridge working 1'30"-2'30", hardening 15-20 min.	5x50 ml
<b>TEMPSI</b>	Silicone-flask separator	60 ml
<b>TDPUN25</b>	Mixing tips	25 pcs

for light curing unit  
LABORLUX3  
see page 58

