## ENA MATRIX | Rubber Dam Instruments Designed by Dr. Lorenzo Vanini

## Iconographic documentation by Dr. Francesco Simoni







EMP5S Ena Matrix Frame for rubber dam small
EMP5M Ena Matrix Frame for rubber dam medium
EMP5L Ena Matrix Frame for rubber dam large





Riubber Dam Punch with 6 holes: the bigger is used for molars on which the clamps are placed; the second for clamp-free molars; the third for premolars and anteriors (the smaller are generally indicated for small teeth and pedodontics).





Forceps for clamps Also suitable also for Ena Matrix Rings



## **RUBBER DAM CLAMPS**







Clamp 22 for superior and inferior molars







Clamp 2 for superior and inferior premolars







Clamps 41 for superior and inferior anteriors, when used in pair ideal for veneers cementation







Clamp 51 small size for superior and inferior anteriors







Clamp 61 medium size for superior and inferior anteriors