




**Product identification:** Smart Propy rubbers and brushes

File N° :0006

<p><b>Manufacturer company:</b> MICERIUM S.p.A.- 16036 Avegno (Ge) Italy – Via G. Marconi, 83 Tel (+39)01857887870 Fax (+39)01857887970 micerium@micerium.it; contabile@pec.micerium.it; www.micerium.it</p>
<p><b>Product reference:</b> ECP01-1, ECP02-1, ECP03-1, ECP04-1, ECP04L-1, ECP31-1, ECP41-1, ECP42-1, ECP43-1 ECP01, ECP02, ECP03, ECP04, ECP04L, ECP31, ECP41, ECP42, ECP43 ECP01B, ECP02B, ECP03B, ECP04B, ECP04LB, ECP31B, ECP41B, ECP42B, ECP43B</p>
<p><b>Description:</b> Medical grade Synthetic Polymer latex free (silicone rubber) and nylon brushes</p>
<p><b>Destination of use:</b> The rotary instruments can be used alone for prophylaxis or in combination with prophylaxis pastes</p>

Micerium Ref.	Description	picture	Material	Size	Colour	Shank	Shape	Stiffness	Sterizable\Disposable
ECP01-1	SMART PROPHY 1pc		nylon brush + Nickel Plated Brass	according to reference sample	pink	RA	CUP	soft	Disposable 134°C min. 3 minutes
ECP02-1	SMART PROPHY 1pc		nylon brush + Nickel plated Brass	according to reference sample	white	RA	CUP	natural medium	Disposable 134°C min. 3 minutes
ECP03-1	SMART PROPHY 1pc		nylon brush + Nickel plated Brass	according to reference sample	blue	RA	CUP	hard	Disposable 134°C min. 3 minutes
ECP04-1	SMART PROPHY 1pc		nylon brush + Nickel plated Brass	according to reference sample: Diam. 2,0 / length 18,5; trim 3,5.	white	RA	Microtuft	Stiff	Disposable 134°C min. 3 minutes
ECP04L-1	"SMART PROPHY RA Small Propy Brush Natural White - 1pc"		nylon brush + Nickel plated Brass	according to reference sample: Diam. 2,0 / length 23,5; trim 3,5.	white	RA	Microtuft	Stiff	Disposable 134°C min. 3 minutes
ECP31-1	PURPLE RUBBER RA 1pc		rubber latex free + Nickel Plated Brass	according to reference sample: Length 22; Cup 10	purple	RA	CUP	medium	Disposable 134°C min. 3 minutes
ECP41-1	SMART PROPHY MINI RA Soft Pink Polisher - 1pc		Medical grade synthetic polymer latex free + Nickel Plated Brass	according to reference sample: Length 18; Cup 6	pink	RA	CUP	soft	Disposable 134°C min. 3 minutes

Micerium Ref.	Description	picture	Material	Size	Colour	Shank	Shape	Stiffness	Sterizable\Disposable
ECP42-1	SMART PROPHY MINI RA Medium Purple Polisher - 1pc		Medical grade synthetic polymer latex free + Nickel Plated Brass	according to reference sample: Length 18; Cup 6	purple	RA	CUP	medium	Disposable 134°C min. 3 minutes
ECP43-1	SMART PROPHY MINI RA Hard Blue Polisher - 1pc		Medical grade synthetic polymer latex free + Nickel Plated Brass	according to reference sample: Length 18; Cup 6	blue	RA	CUP	hard	Disposable 134°C min. 3 minutes

<b>Class:</b> IIa
<b>Composition (main raw material):</b> Synthetic silicone (Medical grade synthetic polymer latex free) or nylon. Metal mandrel JIS H320 C3604 BD or JIS H3250 C3604.
<b>Aspect:</b> Solid
<b>Color:</b> see sheet
<b>Odor:</b> Odorless
<b>Instructions for use:</b> Mount the instrument on a contra-angle handpiece applied to a micromotor. It is advisable to use the product alone or in combination with prophylaxis paste as the last step of the oral hygiene session to smooth the tooth surface after using ablaters. <ul style="list-style-type: none"> <li>• Max speed: 10.000 rpm.</li> <li>• Cooling: use with water spray to reduce heat.</li> <li>• Disposable</li> <li>• Sterilizable at 134°C for 3 minutes before use.</li> </ul>
<b>Warnings and Safety:</b> The product is not dangerous if standard safety measures for rotary dental instruments and for general injury prevention are respected and if the personal protection equipment is used properly (see safety data sheet)
<b>Packaging:</b> Available in packs of 12 pcs . (ECP01, ECP02, ECP03, ECP04, ECP04L, ECP31, ECP41, ECP42, ECP43) and 100 pcs. (ECP01B, ECP02B, ECP03B, ECP04B, ECP04LB, ECP31B, ECP41B, ECP42B, ECP43B). The following codes are used as a single code: ECP01-1, ECP02-1, ECP03-1, ECP04-1, ECP04L-1, ECP31-1, ECP41-1, ECP42-1, ECP43-1
<b>Shelf life:</b> NA
<b>Storage:</b> NA
<b>Min. and Max. Temperature in C°:</b> NA
<b>Interruption at temperature:</b> NA
<b>Transport:</b> NA
<b>Min. and Max. Temperature in C°:</b> NA
<b>Interruption at temperature:</b> NA
<b>Max hours of transport:</b> NA
<b>Freeze resistance:</b> NA

**Warnings for User**

The information included in this document are provided in good faith and based on the present knowledge. They are only indications and in no way binding, in particular for what concerns the violation or the prejudice of other people's rights using our products. Micerium S.p.A. ensures that its products are compliant to their sales specifications. These indications cannot be used as replacement for the necessary preliminary tests that can ensure that a product is suitable for a specific use. The users are responsible to ensure the conformity to the local regulation and to get certifications and authorizations that may be necessary. The users must check that their version of this document is updated, Micerium S.p.A. remains at their disposal for any further information.

Sales Manager Signature:



Data:02/05/2019