



Odt. Daniele Rondoni

Dental Laboratory in Savona | Italy

Modern complex solutions on implants require a special attention not only to the “white” dental tissues but also to the “pink” aesthetic. With the use of “Tender Pink” system I can reproduce with high mimetism the gingival areas: for the different anatomic components I have at my disposal different bodies with calibrated opacity and shades to reproduce the different chromatic combination of the gingival tissues and using the Stains, it is possible to put in evidence also the most fine details.

The use of “Tender Pink” system it is also very interesting in the direct technique by the clinician that can complete directly on the patient that components of structure on implants that cannot be covered in the laboratory, like the screw holes of Toronto Bridge. For the dentist, it is easy to apply the composite Tender Pink and to obtain with simple technique all the mimetisms also for the “pink aesthetic”

Tender Pink

Case Studies Case Studies

