

Ena White



Ena White

System for treatment of pathogen & iatrogenic discoloration

Vital tooth treated with Ena White Power 35% hydrogen peroxide

ENA WHITE POWER

- 35% Hydrogen peroxide for professional treatment in office
- Results in one visit
- Bicomponent syringe activates only the needed quantity for the application
- Neutral Ph 7.3, no enamel damage
- Use with rubber dam or light curing gingival protection Ena Dam



ENA WHITE REGULAR

- 12% Hydrogen peroxide professional treatment in office
- Treatment of non-vital tooth with walking technique
- Safe and effective treatment for sensitive teeth



ENA WHITE LIGHT

- 6% Hydrogen peroxide for cosmetic bleaching
- Home treatment, simple and safe maintenance
- Package with personalized label
- Pleasant freshness: mint flavor!



ENA WHITE

esthetics



Ref	Description	Q.ty
COSM035E	ENA WHITE Power Kit 3.5% Hydrogen peroxide gel Ena Dam paste	5 g 1,5 g
COSM035	ENA WHITE Power Refills 3.5% Hydrogen peroxide gel	2 x 5 g
COSM035T	Ena White automixing tips	50 pcs
COSM012	ENA WHITE REGULAR ENDO Refills 1.2% Hydrogen peroxide gel	2 x 2 ml
COSM015	ENA WHITE LIGHT HOME for Cosmetic Home Bleaching 6% Hydrogen peroxide gel	4 x 3 ml
COSM015B	ENA WHITE LIGHT HOME for Cosmetic Home Bleaching 6% Hydrogen peroxide gel Tray box	2 x 3 ml 1 pc
COSM03B	ENA DAM gingival protection Usable also as block-out	4 g
COSM055	ENA WHITE Desensitizer	2x1,2 ml



DISCOLORATION TREATMENT

ENA WHITE