

ENAMEL PLUS STAIN & INSTRUMENTS


Fluorescent Composite Stains

Unique fluorescent stains used for characterization in both the dental practice and laboratory. Due to high fluorescence, they naturally integrate with the natural tooth. Stains can be used alone or mixed with the composite without changing its fluorescent characteristics.



COSSTAINKIT

MICERIUM BRUSHES FOR COMPOSITE MANIPULATION



Ref	Description	Q.ty
COSSTAINKIT	ENAMEL PLUS STAIN KIT Fluorescent stains in syringe: white, yellow, orange, blue, brown, dark brown Brushes (2 type C - 2 type F - 2 type M)	6x2 g 6 pcs
COSSTAIN1N	White stain refill	2 g
COSSTAIN3N	Blue stain refill	2 g
COSSTAIN4N	Brown stain refill	2 g
COSSTAIN5N	Dark brown stain refill	2 g
COSSTAIN6N	Black stain refill	2 g
COSSTAIN7N	Yellow stain refill	2 g
COSSTAIN8N	Orange stain refill	2 g
COSSTAIN T	Transparent stain refill	2 g
COSSTAINO1	Mixing palette with cover for improving composite modelling	1 pc
COSPEN11	Brush C for stains	1 pc
COSPEN12	Brush F for fossa	1 pc
COSPEN13	Brush M for manipulation	1 pc
TPENS	Soft white silicone brush with small tip for composite handling	5 pcs
TPENN	Hard black silicone brush with small tip for composite handling	5 pcs
TPENMC	Silicone brush blue tip with metal handle 45°/45° sterilizable for composite handling	5 pcs
TPENMS	Silicone brush blue tip with metal handle 90°/180° sterilizable for composite handling	5 pcs
TPENMP	Silicone brush blue tips refill	5 pcs
TLV1	Composite application instrument for anteriors including 1 fixed spatula 1 grey extra flexible spatula 1 white flexible spatula	1 pc 1 pc 1 pc
TLV2	Composite application instrument for posteriors	1 pc
TRAW	White flexible tip	2 pcs
TRAG	Grey ultra flexible tip	2 pcs