TENDER magnetic

TENDER RTRAY Magnetic

PRECISION DIE SYSTEM

This system was born from an idea of a dental technician especially conceived for the dental technician. Thanks to our partnership with mechanical firms and specialised engineers we have been able to create and patent two systems that are compatible and allow to realise any type of dental prosthesis.

Tender Tray "Magnetic" and Tender Flask In Press have been designed with geometrical shapes that are ideal to use with any dental material. The research and the knowledge of the raw materials that are used to make these products are a strong point of our company, whose aim is to offer a unique quality of its kind.

Tender Tray "Magnetic" is a system for sectionable precision models with geometrical and magnetic repositioning. The plate in non-glare material is suitable for Cad-Cam scan.

Tender Flask In Press is an innovative flask with different covers that can hold the Tender Tray "Magnetic", and it can be used with both injection and pressure technique.



Tender Tray 1
with pin retentions and
metal support for "magnetic" dies



Tender Tray 2
with pin retentions
totally made of plastic



Tender Tray 1 and 3 allow the use of a small magnet to be glued under the die to improve the stability and make it incomparable. Moreover the metal support allows the use of another magnet for articulating.



Tender Tray 3
with knurled retentions and metal support for "magnetic" dies on implants



Tender Tray 4 flat without retentions for split cast

Tender Tray 2, is more economical but still has a magnet that facilitates articulating. It can be used as an alternative to the carrying base if the models are sent to dentist on articulator.

Tender Tray 4, ideal for split cast, can be used in the flask same as Tender Tray 1 and 3

TENDER RTRAY Magnetic

FAST, ECONOMICAL AND..... "MAGNETICALLY" ACCURATE

TENDER TRAY Magnetic INTRO KIT

	Tender Tray 1 base Tender Tray 2 base Rubber flap for model casting Rubber box for model casting Aligner / model extractor Magnets for articulating Small magnets for dies Carrying case base	2 pcs 1 pc 1 pc 1 pc 1 pc 2 pcs 10 pcs 2 pcs
	TENDER TRAY 1 <i>Magnetic</i> REFILL Base with pin retentions and metal support	10 pcs
•	TENDER TRAY 2 REFILL Base with pin retentions totally made of plastic	10 pcs
	TENDER TRAY 3 <i>Magnetic</i> REFILL for IMPLANTS Base with knurled retentions and metal support	10 pcs
	TENDER TRAY 4 Magnetic REFILL for SPLIT CAST Base with metal support, no retentions (smooth)	
0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	TX41 Rubber flap for model casting TX42 Rubber box for model casting TX43 Aligner / model extractor TX44 Magnets for articulating TX45 Small magnets for dies TX15 Carrying case bases for Tender Tray 1 / 2 TX35 Carrying case bases for Tender Tray 3	10 pcs 2 pcs 2 pcs 10 pcs 100 pcs 24 pcs 24 pcs




 Ξ

TX10

TX20

TX30

IXF12





TENDER FLASK IN-PRESS

Flask with different covers:

- transparent for composite
- glass fiber for acrylic with silicone mold
- aluminium for acrylic with plaster mold







Toronto bridge
Total prosthesis
Surgical Guide
Temporaries
Bite
Model duplication
Composite cronwn & bridge

TENDER FLASK IN-PRESS

REF. TXF2 SMALL FLASK FOR COMPOSITE

Aluminium base Stainless steel pins "small" 3 pcs Polycarbonate Crystal cover "small" Plates with metal base and magnet 2 pcs Metal disc for magnet 2 pcs Stainless steel latch with click closure Pressure plate Plaster extractors 2 pcs Closing caps for upper holes 6 pcs Filling Spheres 100 pcs



REF. TXF1 COMPLETE FLASK

Aluminium base Stainless steel pins "standard" 3 pcs Polished aluminium cover Glass fiber cover Polycarbonate crystal cover "standard" Stainless steel pins "small" 3 pcs Polycarbonate Crystal cover "small" Plates with metal base and magnet 2 pcs Metal disc for magnet 2 pcs Stainless steel latch with click closure Pressure plate Plaster extractors 2 pcs Closing caps for upper holes 6 pcs Filling Spheres 100 pcs

VERTICOLATORI

Aluminium base with 3 colonnette in acciaio Plates with metal base and magnet Spaziatori varie misure e fermi di chiusura Coperchio verticolatore REF. TX5C VERTICOLATORE x COMPOSITO Verticolatore with Crystal cover REF. TX5B VERTICOLATORE RESINA/CERAMICA



TENDER ~ TRAY Magnetic

STEP BY STEP



Tender Tray 1



Impression alignment



Magnet glued on the dies



Rubber box and flap for model casting



Model removal from rubber box



Sectioned model



Model casting in Tender Tray




Model extrusion



Magnet for articulating



Model casting in the impression



Cutting and finishing



Tender Tray usable in flask

DOUBLE IMPRESSION TECHNIQUE



Aligner for double impression



Trays on articulator using magnets



Positioning of double impression



Models on articulator